

CF-609N

Polyethersulfone Dialysate Filter High Endotoxin Removal Performance

Specifications

- Sterilization : Gamma
- Hollow Fiber
 - Membrane : Hydrophobic Polyethersulfone
 - Membrane Area : 0.6m²
 - Inner Diameter : 500μm
 - Membrane Thickness : 150μm
- Housing : Polycarbonate
- Potting Compound : Polyurethane
- O-Ring : Silicone Rubber
- Maximum Operating Pressure : 1.5kgf / cm²
- High Water Permeability and Pressure Loss (37°C): 0.07kgf / cm² (at 500mL / min)
- Cut-Off Point (Molecular Weight) : 6000
- Package : 12 Pcs / Carton

CE₀₁₂₃

