

Hypodermic Needle

Provides Superior Performance Having Sharp and Smooth Venipuncture

Specifications

- Standard : ISO standard or Nipro standard
- Bevel : Lancet regular or Lancet short bevel
- Wall Thickness : Ultra-Thin Wall(UTW), Thin Wall(TW), or Regular(RW)
- Hub : Transparent color
- Sterilization : ETO / Gamma

CE₀₁₂₃



Product Range

Package : 10,000 Pcs / Carton

Gauge	Color Code	Length (mm)
16	White	25
		38
18	Pink	25
		32
		38
		50
19	Cream	25
		32
		38
20	Yellow	50
		70
		16
		25
21	Deep Green	32
		38
		50

Gauge	Color Code	Length (mm)
22	Black	25
		32
		38
		50
23	Deep Blue	16
		25
		32
24	Medium Purple	38
		25
25	Orange	32
		13
		16
		19
		25
26	Brown	38
		13
		22
27	Medium Grey	13
		19
		32
		38
29	Red	13

* Available upon requests for other sizes, blunt end and Non-sterile.