



# Dental Needle

High Performance with High Quality and Reliability of Products

**Specifications**

- E Type : Thread Screw by mm
- A Type : Thread Screw by inch
- Sterilization : ETO
- Proprietary : Anti-needle breakage design



**Product Range**

Package : 6,000 Pcs / Carton, 10,000 Pcs / Carton

Gauge	Color Code		Needle Length							
	E Type	A Type	(12mm)	(16mm)	(21mm)	(25mm)	(30mm)	(32mm)	(35mm)	(41mm)
25G (0.50mm)	Dark Red	Red			○	○	○	○		
27G (0.40mm)	Yellow	CadmiumYellow			○		○	○	○	○
30G (0.31mm)	Green	Blue	○	○	○	○		○		

\*Other needle length available upon request