



## Blood Tubing Set

### Assuring Compatibility with Most Types of Hemodialysis Machines

#### Specifications

- Bloodline Tube : Medical grade PVC tube (TOTM is available)
- Injection Site : Large finger guard
- Drip Chamber : Available various diameter sizes
- Clamp : Clamps to secure blocking performance
- Recirculation Connector : The simplicity and safety for priming operation
- Remarks : Can be customized with all orders.
- Sterilization : ETO / Steam / Gamma

CE<sub>0123</sub>