

Synthetic Hollow Fiber Dialyzer

- High Flux -

BPA
FREE



CE₀₁₂₃

Our Technology made it Possible

ELISIO™-H

the Next Generation Dialyzer

- ELISIO™-09H
- ELISIO™-11H
- ELISIO™-13H
- ELISIO™-15H
- ELISIO™-17H
- ELISIO™-19H
- ELISIO™-21H
- ELISIO™-25H

PERFORMANCE

* ¹ Clearance(mL/min)	Qb (mL/min)	EL-09H	EL-11H	EL-13H	EL-15H	EL-17H	EL-19H	EL-21H	EL-25H
Urea	200	190	193	196	198	198	199	200	200
	300	246	257	272	278	285	288	291	294
	400	278	298	316	326	337	345	348	363
Creatinine	200	177	184	191	196	197	198	199	200
	300	218	233	250	259	268	273	275	285
	400	242	261	280	296	306	314	326	342
Phosphate	200	163	171	178	184	188	192	195	196
	300	200	213	230	241	254	258	265	276
	400	223	246	265	275	292	305	314	329
Vitamin B ₁₂	200	116	128	140	150	157	164	166	178
	300	134	148	165	180	190	200	206	224
	400	139	161	181	194	211	222	228	247
Inulin	200	81	86	96	102	110	119	124	153
	300	89	94	102	112	121	132	145	171
	400	92	96	109	118	129	139	151	182
Myoglobin	200	58	63	74	84	91	101	104	116
	300	61	68	80	90	98	107	111	126
	400	63	76	84	94	107	113	122	137
* ² KUF (mL/hr/mmHg)		53	59	64	67	74	76	82	93
* ³ KoA Urea		778	916	1190	1351	1614	1771	1976	2269

*¹ In-Vitro Test Conditions(EN 1283): Qd 500mL/min, Qf 10mL/min

*² KUF(EN 1283): Bovine Blood(Hct 32±2%, Protein 60g/L, 37°C), Qb 300mL/min

*³ KoA: Qb 300mL/min, Qd 500mL/min, Qf 10mL/min

SPECIFICATION

	EL-09H	EL-11H	EL-13H	EL-15H	EL-17H	EL-19H	EL-21H	EL-25H	
Effective Surface Area (m ²)	0.9	1.1	1.3	1.5	1.7	1.9	2.1	2.5	
Priming Volume (mL)	62	70	85	95	105	115	130	149	
Effective Length (mm)	212	228	245	259	271	281	290	305	
Inner Diameter (μm)	200								
Membrane Thickness (μm)	40								
Maximum TMP (mmHg)	500								
Pressure Drops	Qb/Qd (mL/min)	200/500							
	Blood/Dialysate (mmHg)	82/19	78/21	76/19	70/21	67/17	67/19	65/17	62/17
Membrane	Polynephron™ (polyethersulfone)								
Housing	Polypropylene								
Header	Polypropylene								
Potting Compound	Polyurethane								
Sterilization Method	Gamma-ray								
Package	24 Pcs/Carton								

NIPRO

Manufacturer site and address:

Nipro Corporation Odate Factory
8-7 Hanukiyachi, Niida, Odate-Shi,
Akita 018-5794, Japan

華江醫療儀器股份有限公司 Hua Chiang Medical Co., Ltd.

[台北公司]
地址:台北市羅斯福路三段245號4樓之2
電話:(02)2362-0396~7
傳真:(02)2362-0398
E-mail:taipei@huachiang.com

[台南公司]
地址:台南市南區新孝路197.199號
電話:(06)264-7878~0
傳真:(06)263-6966
E-mail:tainan@huachiang.com