

Ordering information

Our Continuous Infusion Surefusers™+							
50 mL Surefuser+		100 mL Surefuser+		250 mL Surefuser+			
Code	mL/hr	Code	mL/hr	Code	mL/hr	Infusion time	Colour code
SFS-0505HRP	10,0	SFS-1005HRP	20,0			5 hours	Purple
SFS-0512HRP	4,2	SFS-1012HRP	8,4			12 hours	Pink
SFS-0501DP	2,1	SFS-1001DP	4,2	SFS-10-25P	10,0	1 day	Blue
SFS-0502DP	1,0	SFS-1002DP	2,0	SFS-5-25P	5,0	2 days	Green
SFS-0503DP	0,7	SFS-1003DP	1,4	SFS-3,5-25P	3,5	3 days	Cream
SFS-0505DP	0,4	SFS-1005DP	0,8	SFS-2-25P	2,0	5 days	Yellow
SFS-0501WP	0,3	SFS-1001WP	0,6	SFS-1,5-25P	1,5	1 week	Orange
Pieces per carton: 30		Pieces per carton: 30		Pieces per carton: 15			

Our Short Infusion Time Surefusers™+					
100 mL Surefuser+		250 mL Surefuser+			
Code	mL/hr	Code	mL/hr	Infusion time	Colour code
SFS-1030MP	200,0			1/2 hour	White
SFS-1001HRP	100,0			1 hour	Pale grey
		SFS-250-25P	250,0	1 hour	Light purple
		SFS-125-25P	125,0	2 hours	Salmon pink
		SFS-62,5-25P	62,5	4 hours	Khaki
Pieces per carton: 30		Pieces per carton: 15			

Our Variable Infusion Time Surefusers™+				
100 mL Surefuser+		250 mL Surefuser+		
Code	mL/hr	Code	mL/hr	Colour code
SMA-0123-10P	0,5 - 1,5 - 2,5 - 3,5	SMA-2345-25P	2,0 - 3,0 - 4,0 - 5,0	Green
SMA-2345-10P	2,0 - 3,0 - 4,0 - 5,0	SMA-2468-25P	2,0 - 4,0 - 6,0 - 8,0	
SMA-123-10P	Off - 1,0 - 2,0 - 3,0	SMA-3456-25P	3,0 - 4,0 - 5,0 - 6,0	
Pieces per carton: 30		Pieces per carton: 15		

For information on our products and activities, you can always contact us at the following addresses:



華江醫療儀器股份有限公司
Hua Chiang medical co, ltd

【台北地區】
台北辦事處：台北市羅斯福路三段245號4F-2
電話：(02)2362-0396-7
傳真：(02)2362-0398
Email：taipei@huachiang.com

www.hcmed.com.tw

【台南地區】
台南總公司：702台南市南區新孝路197號
電話：(06)264-7878-0
傳真：(06)263-6966
Email：tainan@huachiang.com



SUREFUSER+™

微量輸液系統



RENAL



Our commitment to strengthen every partnership.

All Nipro renal products are backed by the Nipro Connect Educational Program – offering clinicians product evaluations, in-services, continuing education units (CEU's) and instructional materials.

www.nipro.com

SUREFUSER+之特色及優點

- 輕巧、舒適，便利穿戴
- 經歐盟認證，安全佳、品質佳
- 流量穩定，安全可靠

SUREFUSERTM+

+ 讓醫護人員及病患有更簡易的使用經驗



輕巧、安靜、舒適的特性，
您可輕易的穿戴在身上。


- 協助病患提早離院
- 更佳的舒緩病患疼痛且增加活動度，病患可較早離院
- 無擾人警報聲或提示音
- 藉由無聲的流量控制器確認輸液速度
- 無需裝機保養維修

+ SUREFUSER+可讓每日的使用更為便利

- Surefuser+易於輸液之技術已眾所皆知許久。
- 自供動力、可值得信賴的裝置：
具彈性的加壓(異戊二烯)可維持穩定的流量
- 本產品可輕易裝設及使用，結合便利與效能於一身。



過濾器



+ 經歐盟認證，安全、品質佳

- 堅固的保護球體
- 高品質組件：不含矽膠、不致過敏
- 高效微粒及自我排出的過濾器(0.2 μ)
- 堅固的端口連接器，附加完整的抗扭結(不含DEHP)輸液管
- 平順且設計良好的流量控制器，可確保流量穩定



流量控制器

+ SUREFUSER+改良之優點

- 完整的輸液管
- 輸液殘留量減低至最少
- 減少灌注時間
- 備有背向檢查閥的創新技術
- 產品標籤載明更多產品資訊，
易於辨認(例:批號、有效期…等參考資訊)



+ 可廣泛應用於各式的治療

穩定性資訊：
藥物之穩定性皆由實驗室測試完成，安全可靠。