


ARTERIAL VENOUS FISTULA NEEDLE SET

HIGH QUALITY PERFORMANCE WITH EXTRA THIN WALL NEEDLE



SPECIFICATION


WING : TURNABLE WING OR FIXED WING
 CANNULA : EXTRA THIN WALL WITH OR WITHOUT BACK-EYE
 CONNECTOR & CLAMP :

COLOR	CODE
NEUTRAL	W
RED	R
BLUE	B

TUBE VOLUME :
 300mm LENGTH 2.8ml
 150mm LENGTH 1.5ml

STERILIZATION : ETO Gas
 *Available with Gamma Irradiation

STRUCTURE OF BACK-EYE



PRODUCT RANGE

GAUGE	NEEDLE LENGTH	TUBE LENGTH	WING COLOR
14G	3/4" (20mm) 1" (25mm) 1-1/4" (32mm)	150mm/ 300mm	VIOLET
✓ 15G	✓ 3/4" (20mm) ✓ 1" (25mm) ✓ 1-1/4" (32mm)	150mm/ 300mm	IVORY
✓ 16G	✓ 3/4" (20mm) ✓ 1" (25mm) ✓ 1-1/4" (32mm)	150mm/ 300mm	YELLOW-GREEN
✓ 17G	✓ 3/4" (20mm) ✓ 1" (25mm) ✓ 1-1/4" (32mm)	150mm/ 300mm	ORANGE

CODE	BACK-EYE	CLAMP	WING	
			TURNABLE	FIX
-				○
H	○			○
T			○	
C		○		○
HC	○	○		○
HT	○		○	
TC		○	○	
HTC	○	○	○	


PACKAGE : 500 PCS / CARTON FOR 300mm TUBE LENGTH
 1000 PCS / CARTON FOR 150mm TUBE LENGTH

*SINGLE NEEDLE AVAILABLE WITH NEEDLE LENGTH 1"(25mm)



Nipro disposable arterial venous fistula set

Excellent quality and function from lots of experiences and everlasting studies.




A-V FISTULA NEEDLE


Nipro disposable arterial venous fistula set

●FEATURES

1. Ultra thin wall needle to obtain sufficient blood flow.
2. Safer single-handed connection and disconnection as concavo-convex embossed to avoid slippings in injecting.
3. Prepared a soft tube to absorb shocks to needle point. Leak free construction ensures free from hepatitis.
4. Free from contaminants and infections.
5. Wings are color coded.
6. Luer tapered for easy connection with blood circuit.
7. Back eye type needle available upon request.




●OUT LOOKS



●SPECIFICATION

Gage	Color of Wing	Needle length	Tube length	Capacity
13G	Blue	1"(25mm) or 1 1/4"(32mm)	300mm	2.8mℓ
14G	Purple	" "	" "	" "
15G	Yellow	" "	" "	" "
16G	Green	" "	" "	" "
17G	Red	" "	" "	" "

Siliconized, ultra thin wall needles



NISSHO CORPORATION

3-9-3, Honjo-Nishi, Kita-ku, OSAKA, JAPAN 531

Phone:06-372-2331

Telex:5236888 NISHO J

Cable:"NIPRO MED" OSAKA JAPAN


Facsimile:06-371-7422

Nipro Turnable Wing Arterial Venous Fistula Set




• FEATURES

1. Ultra thin wall needle permitting sufficient blood flow
2. Highly crystallized cannula by newly developed siliconization technique
3. Loaded Mini Clamp for easy and quick handling
4. Luer tapered adaptor for easy connection with blood circuit
5. Back eye type is also available for increasing blood flow.
6. Leak free construction ensures free from hepatitis.
7. Contamination and infection free construction
8. Easily distinguishable by color-coded wings.



• OUTLOOKS



• SPECIFICATION

Item No.	Gauge	Wing Color	Needle length	Standard tube length	Capacity
AVF 13	13G	Blue	1" (25 m/m) or 1 1/4" (32 m/m)	300 m/m	2.2 ml
AVF 14	14G	Purple	"	"	"
AVF 15	15G	Ivory	"	"	"
AVF 16	16G	Green	"	"	"
AVF 17	17G	Orange	"	"	"

• PACKING

ITEM	INNER CARTON	MASTER CARTON
AVF with clamp	50 PCS 230 x 490 x 76 (840 g)	500 PCS 462 x 510 x 400 (10 kg)
AVF without clamp	50 PCS 230 x 490 x 76 (820 g)	500 PCS 462 x 510 x 400 (9.5 kg)

* Standard Export Carton Boxes, durable against water/heat/shock/etc., to be used in delivery.

Remarks

AVF A x B C D E - F G

- | | |
|-----------------------------|-------------------------|
| A: gauge size | E: clamp |
| B: needle length (inch) | C : with clamp |
| C: side hole | nothing : without clamp |
| H : with side hole | F: tubing length (cm) |
| nothing : without side hole | G: connector color |
| D: wing style | W : neutral |
| T : turnable wing | R : red |
| nothing : fixed wing | B : blue |

NISSHO CORPORATION

3-9-3, Honjo-Nishi, Kita-ku, OSAKA, JAPAN 531
 Phone: 06-372-2331
 Telex: 5236888 NISHO J
 Cable: "NIPRO MED" OSAKA JAPAN
 Facsimile: 06-371-7422