



BIOHOLE™ BUTTONHOLE DEVICE

FOR THE BUTTONHOLE TECHNIQUE



BIOHOLE™ BUTTONHOLE DEVICE

FOR THE BUTTONHOLE TECHNIQUE

- Accelerated method to create a buttonhole tunnel track
- Limited use of sharp needles
- Ergonomic design to facilitate grip of the holder
- Single-use
- Not made with latex, DEHP, or PVC
- Available in lengths of 3, 5, or 8 mm

CE 0123
Class IIa

PRODUCT COMPLIANCE

- Complies with the following norms and regulations:
 - EN 556-1
 - EN 556-2
 - EN 1041
 - EN ISO 10993-1
 - EN ISO 10993-3
 - EN ISO 10993-5
 - EN ISO 10993-6
 - EN ISO 10993-10
 - EN ISO 10993-11
 - EN ISO 11137-1
 - EN ISO 11137-2
 - EN ISO 11137-3
 - EN ISO 13485
 - EN ISO 14971
 - EN ISO 15223-1
 - EN 62366-1
 - EN ISO 11737-1

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Japan

STERILIZATION AND SHELF LIFE

Sterilization

Gamma ray



Single-use only



Indicator sticker


"Each outer box (or shipping carton) contains a chemical indicator that ensures sterility. This indicator comes in the form of a yellow sticker that turns red when sterilized."

Shelf life

3 years



MATERIALS USED



Not made with
DEHP, Latex, PVC

PRODUCT RANGE OVERVIEW

Reference code	Length (mm)
BIOHOLE-KIT3	3
BIOHOLE-KIT	5
BIOHOLE-KIT8	8

PACKAGING DETAILS

Reference code	Individual box (2 pieces)				Inner box					Outer box				
	L (mm)	W (mm)	H (mm)	Weight (g)	Quantity	L (mm)	W (mm)	H (mm)	Weight (g)	Quantity	L (cm)	W (cm)	H (cm)	Weight (g)
BIOHOLE-KIT3	10,5	1,7	7,1	8,65	10	18,5	7,6	10,8	143,5	100	42,4	19,9	12,1	1052,5
BIOHOLE-KIT	10,5	1,7	7,1	8,65	10	18,5	7,6	10,8	143,5	100	42,4	19,9	12,1	1052,5
BIOHOLE-KIT8	10,5	1,7	7,1	8,65	10	18,5	7,6	10,8	143,5	100	42,4	19,9	12,1	1052,5

Languages

14 languages on all levels of packaging

English (EN), German (DE), French (FR), Spanish (ES), Italian (IT), Dutch (NL), Bulgarian (BG), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Finnish (FI), Norwegian (NO), Ukrainian (UK)

Transport conditions

Keep away from direct sunlight, excessive humidity, dryness, and vibrations.

Storage conditions

Keep away from direct sunlight, excessive humidity, dryness, and vibrations.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE R	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	✓	✓	✓
	✓	✓	✓

