

NIPRO



BIOHOLE™ Needle System

“百悠活”改良鈕扣式穿刺技術及獨特設計的鈍針



BIOHOLE™ Needle system

減少動靜脈穿刺的疼痛

Biohole Needle System“百悠活”改良鈕扣式穿刺技術是卓越的組合。固定的穿刺通道形成後，您可使用Biohole Needle System，它不像一般AVF銳利，但尖銳程度是足以避免穿刺血管壁較厚的病患時，產生反彈現象。

使用BIOHOLE Needle system“百悠活”改良鈕扣式穿刺技術之優點:

- 減少出血時間
- 使用最小的加壓力便能穿刺，減少病患疼痛
- 延長血管通路之使用壽命
- 減少不當穿刺之現象

BioHole Needle - standard		AVD												
Tubing length (cm) and luer lock		Gauge	14			15			16			17		
		Length (mm)	20	25	32	20	25	32	20	25	32	20	25	32
HTC (back eye)	15R		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	30R		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
TC (Venous)	15R		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	30R		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

BioHole Needle - single			SGD											
Tubing length (cm) and luer lock		Gauge	14			15			16			17		
		Length (mm)	20	25	32	20	25	32	20	25	32	20	25	32
HTC (back eye)	15RB		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	30RB		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

BioHole button	Biohole Kit
----------------	-------------

* Contact your local office for order quantities and availability. All needles are supplied with clamp and coloured turnable wings.