

NIPRO

Nipro Disposable Syringe with needle

Injection System

1cc Insulin
1cc Tuberculin
3cc
5cc
10cc
20cc
30cc
50cc



Safety and accuracy is the concept of NIPRO syringe

Nipro Disposable Syringe with needle

FEATURES

產品特點

- 塑膠注射筒及推桿使用聚丙烯(P.P.)製造，材料安全，美觀透明。
- 推桿停止閥可防止推桿意外掉落的危險。
- 大的手指握柄和拇指接觸區，可使安全緊握，減少在不平衡表面滾動的機會。



SPECIFICATION

capacity needle size		LUER SLIP					LUER LOCK			
		1ML	3ML	5ML	10ML	20ML	3ML	5ML	10ML	20ML
20G	1/2" (13mm)						●			
	1" (25mm)		●	●			●	●	●	
	1-1/4" (32mm)			●			●	●	●	
	1-1/2" (38mm)			●			●	●	●	
21G	1/2" (13mm)				●	●			●	●
	5/8" (16mm)						●			
	1" (25mm)		●	●	●	●	●	●	●	●
	1-1/4" (32mm)			●	●	●	●	●	●	●
	1-1/2" (38mm)		●	●	●	●	●	●	●	●
22G	5/8" (16mm)						●			
	3/4" (19mm)		●				●			
	1" (25mm)		●	●	●	●	●	●	●	●
	1-1/4" (32mm)		●	●	●	●	●	●	●	●
	1-1/2" (38mm)			●	●	●	●	●	●	●
23G	5/8" (16mm)						●			
	3/4" (19mm)		●							
	1" (25mm)		●	●		●	●	●	●	●
	1-1/4" (32mm)		●	●	●	●	●	●	●	●
	1-1/2" (38mm)			●			●	●		
24G	1" (25mm)		●				●	●		
	1-1/4" (32mm)						●	●		
	1-1/2" (38mm)						●	●		
25G	5/8" (16mm)	●	●				●			
	1" (25mm)	●	●	●			●	●		
	1-1/2" (38mm)						●			
26G	3/8" (10mm)	●								
	1/2" (13mm)	●								
27G	1/2" (13mm)	●							●	

NIPRO

Manufacturer site and address:

PT. Nipro Indonesia Jaya
Kawasan Industri Suryacipta, JL. Surya Utama Kav I
22B, 23 & 24, Desa Kutamekar, Kec. Ciampel,
Karawang, Jawa Barat, Indonesia

華江醫療儀器股份有限公司 Hua Chiang Medical Co., Ltd.

[台北公司]

地址:台北市羅斯福路三段245號4樓之2
電話:(02)2362-0396~7
傳真:(02)2362-0398
E-mail:taipei@huachiang.com

[台南公司]

地址:台南市新孝路197號
電話:(06)264-7878~0
傳真:(06)263-6966
E-mail:tainan@huachiang.com