



GLUCOPRO Insulin Syringe

Insulin Syringe with Low Dead Space
Easy to Use and Reliable Quality

Specifications

- Barrel : Polypropylene
- Plunger : Polypropylene
- Gasket : Thermoplastic Elastomer / Rubber Elastomer
- Sterilization : ETO
- Insulin unit : Available to select from 40 units or 100 units



Product Range

Syringe Capacity	Needle		Package
	Gauge	Length	
0.3cc / mL 0.5cc / mL 1.0cc / mL	27G	1 / 2" (13mm)	Wrapper / Blister / Multipack / Pharmacy Pack
		5 / 16" (8mm)	
	28G	1 / 2" (13mm)	
		5 / 16" (8mm)	
	29G	1 / 2" (13mm)	
		5 / 16" (8mm)	
	30G	1 / 2" (13mm)	
		5 / 16" (8mm)	
	31G	1 / 2" (13mm)	
		5 / 16" (8mm)	