

Syringe

Smooth Performance and Transparent with High Grade Material

Specifications

- Barrel : Polypropylene
- Plunger : Polyethylene / Polypropylene
- Gasket : Natural Rubber
- Sterilization : ETO



Product Range

Syringe Capacity	Luer Tip			Individual Packing
	Slip	Lock	Catheter	
1mL U100 Insulin	○			Wrapper or Blister
1mL Tuberculin	○			Wrapper or Blister
3mL	○	○		Wrapper or Blister
5mL	○	○		Wrapper or Blister
10mL	○	○		Wrapper or Blister
20mL	○	○		Wrapper or Blister
30mL	○	○		Blister
50mL	○	○	○	Blister

Product Range (Non-CE marked)

Syringe Capacity	Luer Tip				Individual Packing
	Slip	Lock	Catheter	Eccentric	
1mL	○	○			Wrapper or Blister
3mL	○	○			Wrapper or Blister
5mL	○	○			Wrapper or Blister
10mL	○	○		○	Wrapper or Blister
20mL	○	○		○	Wrapper or Blister
30mL	○	○	○	○	Wrapper or Blister
50mL	○	○	○	○	Wrapper or Blister

Product Range with Needle (Non-CE marked)

O.D Gauge	20G	21G	22G	23G	24G	25G	26G	27G
O.D Metric	0.9mm	0.8mm	0.7mm	0.6mm	0.55mm	0.5mm	0.45mm	0.4mm
ISO Standard Color	Yellow	Deep Green	Black	Deep Blue	Medium Purple	Orange	Brown	Medium Grey
NIPRO Color	Yellow	Deep Green	Black	Light Blue	Red	Blue	Beige	Grey
Length	1/2", 5/8", 3/4", 1", 1-1/4", 1-1/2"							

Needle size mounted on syringe available upon request.