

# SAFETOUCH™ AVF



## Minimizing Needle Stick Injury with Reliable Safety Function

### Specifications

- Wing : Turnable Wing or Fixed Wing
- Cannula : Extra thin wall with/without back-eye

- Tube Length : 270mm  
120mm
- Sterilization : ETO / Gamma

### Connector & Clamp

Color	Code
White	W
Red	R
Blue	B



### Product Range

Gauge	Needle Length	Tube Length	Wing Color		Remarks
			NIPRO	ISO	
14G	1" (25mm)	120(mm) or 270(mm)	Violet	Pale-Green	with Back Eye, Clamp, Turnable or Fixed Wing available upon your request
15G	1" (25mm) 1-1/4" (32mm)	120(mm) or 270(mm)	Ivory / Blue*	Blue-Grey	
16G	1" (25mm) 1-1/4" (32mm)	120(mm) or 270(mm)	Light-Green	White	
17G	1" (25mm) 1-1/4" (32mm)	120(mm) or 270(mm)	Orange	Red-Violet	

Package : 500 Pcs / Carton for 270mm Tube Length. 1000 Pcs / Carton for 120mm Tube Length.

\* for USA market