



NIPRO
EUROPE

SAFETOUCH™ CATH

IV-Catheter with safety device

NEW passive safety technology preventing needle stick injuries

Transparent body for easy flashback confirmation

Secure grips for a precise insertion

Polyurethane catheter, biocompatible and thermo sensitive material,
softens at body temperature

All components are medical grade, latex and PVC free




RECOGNIZED FOR QUALITY, PREFERRED FOR VALUE AND INNOVATION

Product Range

Item number*			Catheter		Needle	Flow Rate	Colour
Straight	Winged	Winged & Port	Size G (mm)	Length Inch (mm)	Size (G)	ml/min	
NIC-SF-1650P-ET	NIC-SF-1650PW-ET	NIC-SF-1650PWP-ET	16G (1.7mm)	2" (50mm)	18G	185 ml/min	Grey
NIC-SF-1832P-ET	NIC-SF-1832PW-ET	NIC-SF-1832PWP-ET	18G (1.3mm)	1¼" (32mm)	20G	93 ml/min	Green
NIC-SF-1850P-ET	NIC-SF-1850PW-ET	NIC-SF-1850PWP-ET		2" (50mm)			
NIC-SF-2032P-ET	NIC-SF-2032PW-ET	NIC-SF-2032PWP-ET	20G (1.1mm)	1¼" (32mm)	22G	56 ml/min	Pink
NIC-SF-2050P-ET	NIC-SF-2050PW-ET	NIC-SF-2050PWP-ET		2" (50mm)			
NIC-SF-2225P-ET	NIC-SF-2225PW-ET	NIC-SF-2225PWP-ET	22G (0.8mm)	1" (25mm)	24G	37 ml/min	Blue
NIC-SF-2232P-ET	NIC-SF-2232PW-ET	NIC-SF-2232PWP-ET		1¼" (32mm)			
NIC-SF-2419P-ET	NIC-SF-2419PW-ET	NIC-SF-2419PWP-ET	24G (0.7mm)	¾" (19mm)	26G	22 ml/min	Yellow
NIC-SF-2425P-ET	NIC-SF-2425PW-ET	NIC-SF-2425PWP-ET		1" (25mm)			

* Individually packed in a blister - inner box 50pcs - export carton 2.000pcs



Nipro Europe NV
European Headquarters
 Weihoek 3h
 1930 Zaventem
 Belgium
 T: +32 (0)2 725 55 33
 F: +32 (0)2 725 70 41
info@nipro-europe.com
www.nipro-europe.com

Nipro Europe - France
 Biôpôle Clermont-Limagne
 63360 Saint Beuzire
 France
 T: +33 473 33 41 00
 F: +33 473 33 41 09

Nipro Europe - Poland
 ul. al. Jana Pawla II 15
 00-828 Warszawa
 Poland
 T: +48 22 653 85 98
 F: +48 22 653 85 06

Nipro Europe - Netherlands
 Regus Business Center
 Verlengde Poolseweg 16
 4818 CL Breda
 Netherlands
 T: +31 (0)76 524 50 99
 F: +31 (0)76 524 46 66

Nipro Europe - Russia
 12 Krasnopresnenskaya
 Nab. Office nr. 1408
 10th Floor - entrance nr. 6
 World Trade Center
 123610 Moscow
 Russia
 T: +7 495 258 1364
 F: +7 495 258 1367

Nipro Europe - Spain
 Poligono Los Frailes n° 93 y 94
 Daganzo
 28814 Madrid
 Spain
 T: +34 91 884 5531
 F: +34 91 878 2840

Nipro Europe - Italy
 Centro Direzionale Milanofiori
 Strada 1 - Palazzo F/1
 20090 Assago - Milano
 Italy
 T: +39 02 57 50 00 57
 F: +39 02 57 51 81 11