

Scalp Vein Set



The Excellent Soft Wing to Skin Surface for Reduced Incompatibility

Specifications

- Cannula : Thin walled stainless needle
- Wing : Color coded flexible wing for easy manipulation
- Tube Capacity : 0.43mL for 300mm tube length
- Connector : Lock hard connector (luer lock)
Soft connector (luer slip) available upon request
- Tube Length : 12" (300mm), available upon request
- Sterilization : ETO



Product Range

Package : 1,000 Pcs / Carton

Wing Color		Needle		Tube Length
NIPRO	ISO	Gauge (mm)	Length	
Pink	Pink	18G (1.2mm)	3 / 4" (19mm)	12" (300mm)
Brown	Cream	19G (1.1mm)		
Yellow	Yellow	20G (0.9mm)		
Green	Deep Green	21G (0.8mm)		
Black	Black	22G (0.7mm)		
Light Blue	Deep Blue	23G (0.6mm)		
Red	Medium Purple	24G (0.55mm)		
Blue	Orange	25G (0.5mm)		
Beige	Brown	26G (0.45mm)		
Grey	Medium Grey	27G (0.4mm)		