



# Safelet Cath

## Intravenous ETFE Catheter

Smooth Puncture with Thin Wall Cannula and Flexible Catheter Allows to Use for Long Period

### Specifications

- Metal needle : Thin walled stainless needle
- Catheter : Flexible ETFE(Ethylene Tetrafluoroethylene) catheter  
Radio-opaque catheter or Non radio-opaque  
Catheter available upon request
- Sterilization : ETO



### Product Range

Package : 1,000 Pcs / Carton

Hub Color		Catheter (ETFE)		Needle
Iso Standard	NIPRO Standard	Gauge (mm)	Length	Gauge (mm)
Orange	Purple	14G (2.2mm)	2" (50mm)	16G (1.6mm)
Grey	Yellow Green	16G (1.7mm)	2" (50mm)	18G (1.2mm)
Grey	Yellow Green	16G (1.7mm)	2-1/2" (63mm)	18G (1.2mm)
Green	Pink	18G (1.3mm)	2" (50mm)	20G (0.9mm)
Green	Pink	18G (1.3mm)	1-1/4" (32mm)	20G (0.9mm)
Pink	Yellow	20G (1.1mm)	2" (50mm)	22G (0.7mm)
Pink	Yellow	20G (1.1mm)	1-1/4" (32mm)	22G (0.7mm)
Blue	Black	22G (0.9mm)	1-1/4" (32mm)	24G (0.55mm)
Blue	Black	22G (0.9mm)	1" (25mm)	24G (0.55mm)
Yellow	Red	24G (0.7mm)	1" (25mm)	26G (0.45mm)
Yellow	Red	24G (0.7mm)	3/4" (19mm)	26G (0.45mm)