



SUREFUSER™ +

Elastometric Disposable Infusion pump

Excellent Flow Rate for each Therapy and Lightweight Disposable Design

Specifications

- SUREFUSER + is indicated to use for infusion of medications at a continuous and accurate predetermined flow rate. The most important change we have implemented with the SUREFUSER+ is the integration of the infusion line with the balloon. With the new feature, it is unnecessary to attach the infusion line while filling the medication and during the filling operation.
- Sterilization : ETO



Package : 50 / 100 mL Type : 30 pcs / Carton, 250 mL Type : 15 pcs / Carton

Product Range

| Model | Volume | Infusion Rate | Infusion Time | Color Code (Casing) |
|--------------|------------------------------|----------------------|----------------------------|---------------------|
| SFS-0505HRP | 50mL | 10.0 mL / hr | 5 hours | Purple |
| SFS-0512HRP | | 4.2 mL / hr | 12 hours | Pink |
| SFS-0501DP | | 2.1 mL / hr | 1 day | Blue |
| SFS-0502DP | | 1.0 mL / hr | 2 days | Green |
| SFS-0503DP | | 0.7 mL / hr | 3 days | Cream |
| SFS-0505DP | | 0.4 mL / hr | 5 days | Yellow |
| SFS-0501WP | | 0.3 mL / hr | 1 week | Orange |
| SFS-1030MP | | 100mL | 200 mL / hr | 1/2 hours |
| SFS-1005HRP | 20.0 mL / hr | | 5 hours | Purple |
| SFS-1012HRP | 8.3 mL / hr | | 12 hours | Pink |
| SFS-1001DP | 4.2 mL / hr | | 1 day | Blue |
| SFS-1001HRP | 100 mL / hr | | 1 hour | Pale Grey |
| SFS-1002DP | 2.1 mL / hr | | 2 days | Green |
| SFS-1003DP | 1.4 mL / hr | | 3 days | Cream |
| SFS-1005DP | 0.8 mL / hr | | 5 days | Yellow |
| SFS-1001WP | 0.6 mL / hr | | 1 week | Orange |
| SFS-250-25P | 250mL | | 250 mL / hr | 1 hour |
| SFS-125-25P | | 125 mL / hr | 2 hours | Salmon Pink |
| SFS-62.5-25P | | 62.5 mL / hr | 4 hours | Khaki |
| SFS-10-25P | | 10 mL / hr | 1 day | Blue |
| SFS-5-25P | | 5 mL / hr | 2 days | Green |
| SFS-3.5-25P | | 3.5 mL / hr | 3 days | Cream |
| SFS-2-25P | | 2 mL / hr | 5 days | Yellow |
| SFS-1.5-25P | | 1.5 mL / hr | 1 week | Orange |
| SMA-0123-10P | | 100mL | 0.5, 1.5, 2.5, 3.5 mL / hr | — |
| SMA-2345-10P | with Variable Flow Regulator | 2, 3, 4, 5 mL / hr | — | Green |
| SMA-123-10P | | Off, 1, 2, 3 mL / hr | — | Green |
| SMA-2345-25P | | 250mL | 2, 3, 4, 5 mL / hr | — |
| SMA-2468-25P | with Variable Flow Regulator | 2, 4, 6, 8 mL / hr | — | Green |
| SMA-3456-25P | | 3, 4, 5, 6 mL / hr | — | Green |