



# SAFETOUCH™ SCALP VEIN SET

## With Safety Device

Nipro SafeTouch™ Blood Collection Set is specifically designed to prevent needle stick injuries.

Needle tip protection is achieved by a simple procedure which ensures that the needle can never be re-exposed.

- **One hand operation**

Activation of safety mechanism is achieved using one hand only.  
Use the other hand to apply pressure to venipuncture site, as usual.

- **No exposure of sharp needle**

Safety mechanism is activated while needle is being withdrawn. As it is not necessary to withdraw needle from patient, injury risk is eliminated.

- **Secure and user friendly safety mechanism**

Before activation, the mechanism stays securely locked. SafeTouch™ is a user friendly safety mechanism and requires no change of technique. After activation, SafeTouch™ ensures that it cannot be activated inadvertently.

- **Large wings for easy and firm grip**


Large wings are designed to be folded easily and gripped firmly with interlink connection pads. This helps avoid accidental needle stick prior to venipuncture.



RECOGNIZED FOR QUALITY, PREFERRED FOR VALUE AND INNOVATION

### SAFETY MECHANISM

1. Push



2. Pull



3. Click



### PRODUCT RANGE

	Product	Gauge	Needle length	Tubing	Hub color
Luer Slim (Soft connector)	18Gx3/4"	18G	19mm	30cm	Pink
	19Gx3/4"	19G	19mm	30cm	Cream
	20Gx3/4"	20G	19mm	30cm	Yellow
	21Gx3/4"	21G	19mm	30cm	Green
	22Gx3/4"	22G	19mm	30cm	Black
	23Gx3/4"	23G	19mm	30cm	Blue
	24Gx3/4"	24G	19mm	30cm	Violet
	25Gx3/4"	25G	19mm	30cm	Orange
	26Gx3/4"	26G	19mm	30cm	Brown
27Gx3/4"	27G	19mm	30cm	Gray	
Luer Cap (Hard connector)	18Gx3/4"	18G	19mm	30cm	Pink
	19Gx3/4"	19G	19mm	30cm	Cream
	20Gx3/4"	20G	19mm	30cm	Yellow
	21Gx3/4"	21G	19mm	30cm	Green
	22Gx3/4"	22G	19mm	30cm	Black
	23Gx3/4"	23G	19mm	30cm	Blue
	24Gx3/4"	24G	19mm	30cm	Violet
	25Gx3/4"	25G	19mm	30cm	Orange
	26Gx3/4"	26G	19mm	30cm	Brown
27Gx3/4"	27G	19mm	30cm	Gray	

### PACKAGE

Package	Quantity
Inner box	50 pcs.
Outer carton	500 pcs.

Authorized Representative :  
**NIPRO EUROPE NV**  
 Weihoek 3 H, 1930 Zaventem, Belgium  
 Phone: +32 (0) 2 725 55 33 · Fax: +32 (0) 2 725 70 41

Manufacturing facility :  
**Nipro (Thailand) Corporation Limited**  
 10/2 Moo 8, Bangnomko, Sena, Phra Nakhon Si,  
 Ayutthaya 13110, Thailand